

Büchi «flexyclave»

The automated hydrogenation solution for laboratory and scale-up



- Easy handling
- Safe operation
- High reproducibility
- Increased productivity
- Recipe control
- Flexibility
- Graphical data evaluation
- Automatic protocol
- Compact design
- Turn-Key solution

in cooperation with
SYSTAG

BÜCHI – THE WAY TO GET RESULTS!


«flexyclave»

The Turn-Key Hydrogenation Solution Result of a Büchi – Systag cooperation

The easy to use recipe control software with predefined steps and operating procedures makes hydrogenations accurate, reproducible and safe. Processes can be standardized which facilitates up-scaling from laboratory to pilot plant and production scale. Easy operation allows the user to carry out hydrogenations quickly and safely. The manual mode guarantees high flexibility at

all times. Sophisticated alarm- and supervisory functions, allows the system to be run unattended. All process parameters, operator manipulations and results are continuously and automatically recorded, ideal for process evaluation and fine tuning for a next run.

«ecoclave» pressure reactor with glass vessel and gas dosing unit




Visual observation of reaction with glass vessel



Polycarbonate protection door for optimal user safety



Turbine stirrer with gassing stirrer shaft for effective gas dispersion






Including piping and wiring



Compact gas dosing unit containing electronics, flow measuring, control and dosing valves

FlexySys-Software

Automation with recipes – Simple operation



Intuitive plant synoptic, because ...
... not only specialists do hydrogenations


Predefined operating procedures for simple and safe operation

- reaction temperature
- reaction pressure
- end criteria

... ready, steady, go


Clear illustration of system status, then ...

... you know where your process is – exactly!



Functions

- ✓ Hydrogenation with MFC
- ✓ Different end criterias
- ✓ Temperature control
- ✓ Reactor / Jacket / Ramping
- ✓ Pressure- / Vacuum control
- ✓ Automatic inerting
- ✓ Automatic aeration
- ✓ Stirrer speed control
- ✓ Recipe control
- ✓ Visual process control using online graphs
- ✓ Data report
- ✓ Offline-calculations



Time	Event
40.007	Phase 1 canceled (Pressure)
40.007	Phase 2 loaded (Saturation)
40.007	Phase 3 loaded (Vent)

Pres End value = 1.200
 Pres Ramp time = 2.0min
 Pres Regulation state = läuft
 Pres Ramp gradient = 19.40bar
 Pres Regulation state = läuft, bee
 Phase 1 canceled (Pressure)
 Pres Regulation state = gestoppt
 Phase 2 loaded (Saturation)
 Timer delay = 1min
 Timer state = finished
 Phase 3 loaded (Vent)
 End value = 5.000

Parameter	Value
Temperature safety	200.0 °C
Pressure safety	30 bar
Stability band	6.3 bar
Master P band	234.5 bar
Slave D band	6.4 bar
Slave P band	8.2 bar
Slave I part	32.2 s

End value = 5.000
 Master P band = 234.5 bar
 Master I part = 50.5 s
 Slave D band = 6.4 bar
 Slave P band = 8.2 bar
 Slave I part = 32.2 s
 Pressure parameters
 Stability band = 6.3 bar
 Stability time = 1.0 min
 P band = 6 bar
 I part = 400.0 s
 D part = 0.0 s
 Factor = 1.00

Editable online trend, then ...
 ... one picture says more than a thousand words

Automatic and continuous documentation, then ...
 ... what is not written down didn't happen

System parameters at a glance, then ...
 ... safety is when you know what's going on

Technical data

PC-Hardware / Software

CPU	®Core™ 2 Duo E7500 with VT (2.93GHz, 3M, 1066MHz)
RAM	min. 2048MB
LAN	2 x Ethernet (100 MBit / s)
Monitor	min. TFT 19" (resolution min. 1280 x 1024)
Operating system	Windows XP Professional
.Net Framework	Microsoft .Net Framework 2.0
Office	Office 2003 or 2007 Basic (Word & Excel)

Reactor

Frame	stainless steel construction with polycarbonate-safety screen
Agitator	«cyclone 075» magnetic drive, 100...3000 rpm, max. torque 75Ncm, gassing stirrer shaft with turbine stirrer
Cover plate	with 7 openings, 5x ¼" NPT, 1x ⅜" NPT, 1x ½" NPT
Vessel	glass vessel Type 1B, 1.0 liter, with jacket and Torion bottom valve, max. -1...+6 bar, max. -50...+200°C
Accessories	Manometer, bursting disc, pressure transducer, Pt100 sensor, speed sensor, 2 control valves, addition funnel with screw cap, baffle
Material	all product touched parts are made of borosilicate glass 3.3 or stainless steel 1.4435 / 1.4571 (316 / 316Ti)
Design	manufactured and tested according PED, AD2000
Dimensions	H x W x D 1000 x 500 x 500 mm

Thermostat is not included. We recommend models from Huber or Julabo. Others upon request

Options

Special executions are no problem. Upon request we can accommodate almost every wish:

- higher pressure and temperature
- other materials (Hastelloy, titanium, etc.)
- bigger reactors (up to 20 liter)
- agitator cyclone 300
- liquid dosage
- pH measurement / control
- other applications
- CFR21 Part 11
- qualification
- PAT – RTA
- custom software functions

Ask us!

Control unit

PLC	2 x 65HC11-Processors, 16 bit A/D converter, RS-232 interface
Operating temperature	10 .. 35 °C
Mains supply	230 VAC, 10A, 50 / 60 Hz, single phase, uninterrupted
Dimensions	H x W x D 510 x 510 x 280 mm
Utilities	Gas supply on top, signals and power on left or right
Online recording	of temperature setpoint, reactor internal temperature, reactor jacket temperature, stirrer speed, hydrogen flow, hydrogen consumption
Safety supervision	of temperature, pressure and hydrogen leakage inside control unit
Continuous recording	of all manual and recipe interventions
Temperature meas. range	-100°C to 300°C, resolution 0.1K
Hydrogen	minimum flow 0.16 to 8 ml per minute maximum flow 0.16 to 25 L per minute maximum pressure 64 bar

